

Pharma  Plus



Microscopy Based Particle Size Analyser

www.promptplus.in



Microscopy Based Particle Size Analyser

Microscopy Based Particle Size Analyser


The microscope has always been the judge in particle characterization. Since it is accepted as the most direct measurement of particle size.

We have automated manual microscopy to replace a tedious, somewhat subjective measurement with a sophisticated technique for quantifying size of a sufficient number of particles to assure statistical confidence with the end result.

Pharma i Plus.....

This tool gives the capability to analyse the particles in the range of 0.5 μm to 2000 μm . The system provides high resolution, accurate and reproducible results. It can measure a wide range of particles and can detect small differences in samples. The new software provides the data you need with only a few clicks.

Pharma i Plus can analyse Dry, Aqueous, Non – Aqueous, Emulsions, Aerosols, Actuators Orifice, Capsule Soft Gelatin, Capsule thickness, Tablet Coating and many more...








The Systems Features

Three Reporting Options in One software :

- i. Automated Particle size counting
- ii. Morphometry (Manual Measurement Length & Area)
- iii. Image Documentation



Automated Particle size counting

Morphometry

Image Documentation



Automated Particle Size Counting


- Advanced processing algorithms for detecting different type of particles.
- Standard and user definable method files to handle different sample types. Easy Creation of product wise Macro.
- Classification of particles on the basis of Major Axis.
- Can process a batch of greater than 50 images at a time for automatic multiple field analysis.
- Report generated in the form of PDF file.
- In built report generator for generating custom reports complete with Data tables, histograms and images.
- Each Particle's data can be transferred to PDF.
- All data stamped with date, time, magnification and calibration information.
- Results available in form of D10, D50, D90, D100 formats or customized as per user requirement.
- Percentage wise classification for Particle size distribution
- Particle with or without agglomerate are separately defined in the result table.
- Report with total Particle count, Min. size, Max. size and Average size of particle

Sample Detection :

Simple and easy methods are available on the go to analyse different samples. Pharma i Plus system has an automated technique which refrains the user from venturing into tedious thresholding techniques.



Captured Image



Detected Image



Reports

Classification (Single Particle) :

User defined classification can be done below :

Class	Count	%Count
0 to 10	1654	65.582871
10 to 25	668	26.486915
25 to 50	182	7.216495
50 to 100	16	0.634417
100 Above	2	0.079302

Results :

Available in D10, D50, D90, D100 formats or customized

Report	Count
D10	3.736921
D50	7.473841
D90	22.421524
D100	106.513694

Result Summary :

Min. Size, Max. Size and Average Size details Unit, Total Particle Count and Field details

Result Summary	(Major Axis)
Unit	Um
Min	1.67
Max	170.9
Average Size	17.02
Field	2
Total Count (with agglomerate)	2560
Total Count (particle)	95


Classification (with Agglomerate):

A separate data for agglomerate as below :

Class	Count	%Count
0 to 10	1805	54.980201
10 to 25	1008	30.703625
25 to 50	378	11.513859
50 to 100	85	2.589095
100 Above	7	0.21322

Graphical Data :


The data is also represented in graphical format





Morphometry

This advanced Software allows you to Measure Particle, Objects Manually & Generate Report. Two Options available in software for Manual measurements. that is Length & Polygon.



Captured Image



Measured Image

Measurements :

Length Option of Length allows to make measurements of Particles/Objects by dragging cursor on particles edge.

Measurements : Length	
Sr.	Length
1	14.19
2	12.78
3	16.08
4	13.29
5	16.83
6	5.39
7	15.19
8	11.17
9	13.91
10	2.02

Result Summary :

Min. Size, Max. Size and Average Size details, Unit, Total Measurements

Result Summary : (Unit – um)	
Measurement	Length
Total Measurement	30 Nos.
Minimum Size	2.02
Maximum Size	20.45
Average Size	7.65

Image Doc

Software allows you to Capture Microscopy Images & Generate Report with individual comments for Image.

Image Doc option useful to capture and record images of Crystalline Sample, Emulsions, Microbiological Slides, Biological Slides, Petri Dish Images.

Captured Image 1



Image 1

crystalline image under polarizing filter

Captured Image 2



Image 2

Particle Image under Brightfield microscope



Microscopy Based Particle Size Analyser



Software Technology

- Developed with Current Modern Technology
- Microsoft SQL
- PDF
- Framework 4.5

System Requirements

- Windows 7, 10 OS with 64 bit
- Minimum : RAM 8 GB, Graphics 2GB, Hard Disk 500GB
- Active USB port 2nos.
- Microscope with Built_In camera with USB output

Key Functions

- Sample ID creation & multiple analysis under same sample id
- Three reporting options in one software
 - i. Automated particle size counting
 - ii. Morphometry manual measurement
 - iii. Image documentation
- Automated particle size counting in D10, D50, D90& D100 format or as per your requirement with classification table, min./max./average size, with agglomerate & without agglomerate, etc
- Automated particle size counting of single and multiple images
- Customized macros
- Manual measurements of length and polygon
- Option to generate report for microbiological & crystallinity images with comments & remarks
- Easy Calibration & alert for calibration expired
- Usage Log with analysis start/end, reviewed by, date & time information

21 CFR PART 11 COMPLIANCE

- User management
- Creation of groups (customized) i.e. Admin, analyst, reviewer, etc
- User creation & group wise rights allocation
- Login policies, Reset password, Change password
- Audit trails with date & time, User details, Sample ID & activities performed with remarks
- Sample ID creation & multiple analysis under same sample ID
- Non-editable report with E-signature & comments of analyst & reviewer
- Printing of reports, Audit trial, List of users, Group wise rights & Usage log with stamp of date & time, Printed by, Software name & version
- Facility to set Data path & Restoration of data



Microscopy Based Particle Size Analyser

Developed by :

Pulse Life Science

B/107, Vrindavan B.Morwala CHS,
Majaswadi, Near Gajanan Temple,
Jogeshwari (E), Mumbai -400060.
Mob.: 9324479648/9022592910
Email : pulselifescience@gmail.com

Marketed by :

Prompt Plus

2, Milap Apt., Navagaon,
Dahisar (w), Mumbai – 400068.
Mob.: 9869073507 /9320073507
Email : prompt_plus@hotmail.com
Website : www.promptplus.in