

## INTENDED USE

RNA WRAPR Cryovial & RNA WRAPR multi collect media provides the ideal solution for long-term sample preservation for large scale biorepositories.

## SUMMARY AND EXPLANATION

RNA WRAPR Cryovial & RNA WRAPR multi collect media provides the ideal solution for long-term sample preservation for large scale biorepositories. The tube and media guarantees the safe storage of your high-value genomic material and samples for future downstream processing without compromising the quality of genetic material.

## PRODUCT FEATURES

- Accurate & reproducible DNA/RNA yields for down-stream processing and analysis.
- Indefinite storage at -20°C down to vapour phase of liquid nitrogen temperatures.
- Free of DNase & RNase.
- Ultra-pure, medical grade polypropylene tubes.
- Easy sample processing.
- Zero contamination shelf-life for up to 2 years.

## PRECAUTIONS

1. Do not use after expiration date.
2. Product is intended for use as supplied. Do not dilute or add other components to RNA WRAPR Cryovial & RNA WRAPR multi collect media.
3. Overfilling or underfilling of tubes will result in an incorrect sample to stabilization solution ratio and may lead to incorrect analytic results or poor product performance.

## CAUTION

- During operation, always wear a lab coat, disposable gloves, protective goggles and mask.
- All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection.
- Products should be disposed of with infectious medical waste.

## RNA WRAPR CRYOVIAL VOLUME RECOMMENDATIONS

Sample Type	Weight/Vol of Sample	Vol of RNA WRAPR
Bacterial Cells	10 <sup>8</sup> cells	1 ML
Animal Cells	10 <sup>7</sup> cells	0.5 ML
Animal tissue	>200 mg	1 ML
Plant tissue	< 500 mg	1 ML
Environment samples	< 200 mg	1 ML

## STORAGE AND STABILITY

Temperature	Sample Stability
-20°C and below	Indefinite

- When stored at room temperature empty RNA WRAPR Cryovial is stable through expiration date.

## LIMITATIONS

1. For single use only.
2. Samples drawn in other preservatives may cause precipitation in RNA WRAPR Cryovial & RNA WRAPR multi collect media.

## ORDERING INFORMATION

Please call us at +91 8088747968 or mail at [hello@azooka.life](mailto:hello@azooka.life) for any queries or assistance.

Additional information can be found online at [www.azooka.life](http://www.azooka.life)

**Cat No. RW103**

**MANUFACTURED AT:**

# 1A, Kushal Garden Arcade, 'C' Block, 5th Floor, Peenya Industrial Area, 2nd Phase, Bengaluru, Karnataka, India- 560058

**ISO 13485 CERTIFIED**

**CE IVDR**