

INTENDED USE

The mWRAPR Cervical Sample Collection Tube is a user-friendly cervical collection kit to conveniently collect, preserve, and transport cervical specimens containing nucleic acids.

SUMMARY AND EXPLANATION

mWRAPR Cervical Sample Collection Tube ensures sample stability during transport/storage at ambient temperature. mWRAPR Cervical Sample collection tube's combination of preservatives and metabolic inhibitors potentially inactivates infectious agents and nucleases, thereby preserving the genetic material for downstream applications (molecular diagnostics, genomics etc). This collection kit can be used for stabilization of DNA and RNA for microbiome profile analysis.

Samples stored in these mWRAPR Cervical Sample Collection Tubes are stable at ambient temperature, and can be frozen for longer-term storage.

PRODUCT FEATURES

- Accurate & reproducible DNA/RNA yields for down-stream processing and analysis.
- Stabilizes DNA for up to 30 days at room temperature
- No cold storage required
- Easy sample processing
- Zero contamination shelf-life for up to 2 years

PRECAUTIONS

- 1.Do not use tubes after expiration date.
- 2.Product is intended for use as supplied. Do not dilute or add other components to mWRAPR Cervical Sample Collection tubes.
- 3.Overfilling or underfilling of tubes will result in an incorrect sample-to-stabilization solution ratio and may lead to incorrect analytic results or poor product performance.


CAUTION

- During operation, always wear a lab coat, disposable gloves, protective goggles and mask.
- All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection.
- Products should be disposed of with infectious medical waste.

DIRECTIONS FOR USE

- To collect cervical samples, insert the speculum into the vagina and locate the cervix at the end of the endocervical canal.
- Insert the cervical swab into the endocervical canal and rotate the swab (360 degrees) for 10-20 seconds to collect the sample.
- Then, withdraw the swab carefully without touching the skin outside the vagina and place the swab into the collection tube containing the preservative solution.
- Break the swab tip at the score line and close the collection tube tightly.

- Once this is done, ensure that the collection tube is tightly closed and invert the tube gently a few times.



STORAGE AND STABILITY

Temperature	Sample Stability
-20°C and below	Indefinite
4°C	Up to 1 year
Ambient Temperature	Up to 30 days
35-40°C	Up to 14 days
50-65°C	Up to 7 days

- When stored at room temperature, empty mWRAPR Cervical Sample Collection tube is stable through expiration date.

LIMITATIONS

1. For single use only.
2. Samples added in other preservatives may cause precipitation in mWRAPR Cervical Sample Collection tubes.

ORDERING INFORMATION

Please call us at +91 8088747968 or mail at hello@azooka.life for any queries or assistance.
Additional information can be found online at www.azooka.life

Cat no. MW-AZ0026
MANUFACTURED AT:

1A, Kushal Garden Arcade, 'C' Block, 5th Floor, Peenya Industrial Area, 2nd Phase, Bengaluru, Karnataka, India- 560058

ISO 13485 CERTIFIED

CE IVDR