

INTENDED USE

The mWRAPR Skin Microbiome Collection Tube is a user-friendly skin microbiome collection kit for stabilizing the nucleic acid content of skin microbiome.

SUMMARY AND EXPLANATION

Conveniently collect, store and transport skin microbiome at ambient temperature with mWRAPR Skin Microbiome Collection Tube for efficient analysis.

mWRAPR Skin Microbiome Collection Tube effectively inactivates pathogens and arrests their growth while preserving pathogenic genetic material for a variety of tests.

PRODUCT FEATURES

- Accurate & reproducible DNA/RNA yields for down-stream processing and analysis.
- Stabilizes DNA for up to 30 days at ambient temperature.
- No cold storage required.
- Easy sample processing.
- Zero contamination shelf-life for up to 2 years.

PRECAUTIONS

1. Do not use tubes after expiration date.
2. Product is intended for use as supplied. Do not dilute or add other components to mWRAPR Skin Microbiome Collection Tube.
3. Overfilling or underfilling of tubes will result in an incorrect sample-to-stabilization solution ratio and may lead to incorrect analytic results or poor product performance.

CAUTION

- During operation, always wear a lab coat, disposable gloves, protective goggles and mask.
- All biological specimens and materials coming in contact with them are considered biohazards and should be treated as if capable of transmitting infection.
- Products should be disposed of with infectious medical waste.


DIRECTIONS FOR USE

Pre-Collection Instructions

- Do not wash your face before sample collection, so the patch captures sample when your skin is in its most natural state.


Swab-based

- Moisten the swab with SCF-1 Solution.
- Hold the swab parallel to the skin surface and firmly rub it back and forth 20-30 times while continuously rotating the swab head.
- Insert the swab head into the mWRAPR stabilization solution tube.
- Break the swab tip at the score line and close the collection tube tightly.
- Vortex well before the isolation.



Patch-based

- Simply stick the Skin patch on your forehead (adhesive side down) for 10 seconds.
- Slowly plunge the skin patch into the mWRAPR stabilization solution in the Tube. Turn over the capped tube and quickly shake it 10 times.
- Vortex well before the isolation.



STORAGE AND STABILITY

Temperature	Sample Stability
-20°C and below	Indefinite
4°C	Up to 1 year
Ambient Temperature	Up to 30 days
35-40°C	Up to 14 days
50-65°C	Up to 7 days

- When stored at room temperature, empty mWRAPR Skin Microbiome Collection Tube is stable through expiration date.

LIMITATIONS

1. For single use only.
2. Samples added in other preservatives may cause precipitation in mWRAPR Skin Microbiome Collection Tube.

ORDERING INFORMATION

Please call us at +91 8088747968 or mail at hello@azooka.life for any queries or assistance. Additional information can be found online at www.azooka.life

Cat No. MW-AZ0022

MANUFACTURED AT:

1A, Kushal Garden Arcade, 'C' Block, 5th Floor, Peenya Industrial Area, 2nd Phase, Bengaluru, Karnataka, India- 560058

ISO 13485 CERTIFIED

CE IVDR