Meditec 1700 Ventilator





Ideal for use in an ICU, Intra-Hospital transport, Air and Ground transport. Simplicity, In-built air supply, versatility and affordable cost



The Meditec 1700 is designed to work in a highly versatile environment, whether you need Invasive or Non-Invasive ventilation, it delivers both. Meditec 1700 is micro-processor controlled design with wide range of Ventilation modes to suit individual preferences.



Features

- Suitable for adult and paediatric patients
- Built-in turbine technology
- Pressure & Flow waveforms with Trends Facility
- FiO₂ Measurement
- Suitable with both dual as well as single limb breathing circuits
- Patient History storage capability
- Automatic Parameters adjustment of patient according to Height & Weight
- Touch screen as well as knob operated
- Provision for 100% oxygen for a period of 2 min before disconnecting for suctioning
- Inspiration hold: Measures patient's lung compliance and resistance, Elasticity, Time constant
- Key lock: Key & LCD Radiogram Touch Screen Lock
- Waveform freeze: Freezes the waveform
- Alarm silence: Mutes audible alarms for 2 minutes
- Alarm reset: Clears visual indicators and messages
- Easy detaches from trolley and moves with patients
- Mobile stand with quick-release
- APNEA backup with option PACV, VACV with pressure limit
- Shock Resistant (Drop Test compliant)
- Altitude compensation





Waveforms



Event Log



Graphics - Wavegraph, Trend & Loop

			11.7	(1:19) 3.3E	READY	*	Manual Inspiratior
• Monit	oring		Table	1 Table2	2 X	pinsp 80 5	отна 16
	Paw			FPE	AK	RATE	EPN
	P PEAK		VTe				- 15
	P MEAN		VE MIN			PEEP	cmH20
	P PAUSE		VE min SPONT				3
	PEEPH	9	RATE	9 02			- u
	CmHbO PEEP		R SPONT	RA		02	54
	CmH2O		BPM Depi	omH			- 21
0.0	CmH2O		biminimi	mile			~ '
	P01 CmHzO		EXEND FLOW	WOE	3v	Ala	arms
Monitoring	Grapi	nics	Tools	Events	System	◄ 3)	i 🔒

Monitoring Parameters





Non-Invasive Ventilation

Meditec 1700 comes in a full range of clinical requirements with invasive ventilation, automated ventilation and non-invasive mask ventilation.

Performance Data

- Valve response time $\leq 5 \text{ ms}$
- Control principle Time cycled, volume constant, pressure-controlled
- Safety relief valve 100 mbar/cmH₂O
- Automatic Leakage compensation
- Breathing hose & Dynamic compliance compensation
- Max. flow for pressure support and spontaneous breathing 180 L/min

Operating Data

- Power consumption Approx. 150VA
- O₂ supply pressure range : 2.4 6.9 bar / 35-90 psig

Weaning

• O₂ Streaming assisted therapy for weaning patients

Invasive Ventilation

- Exsens & Ex-trig value, W.O.B. of patient
- SIGH with interval 30, 60, 90, 120 / breath
- BTPS reduces Tidal Volume error

The 1700 ventilator combines reliability and high performance.

12.1" TFT Graphical Display Unit (Optional)



Monitoring Parameters		
$V_{\rm E}$ Tidal	P _{peak}	
V _E Min	P _{mean}	
Resp	V, Tidal	
P _{Pause}	F _{peak}	
PEEP	O ₂ %	

Optional Parameters			
PR	EtCO ₂	Resp	
SpO ₂	iCO ₂		
Waveforr	n		
P-T	F-T	V-T	

Technical Specifications

_ |

Display Screen	8" Screen with Touch Control & Rotary Knob
Parameters	Set Parameters / Measured Parameters / Alarm Status / I:E Ratio
Waveforms	Pressure-Time, Flow-Time, Volume-Time
Loops	Pressure-Volume, Flow-Volume, Pressure-Flow
Modes	PACV, PSIMV, VACV, VSIMV, AUTO (Emergency Quick Start), Spont+, PRVC, ACPAP+, PRVC+SIMV
Tidal Volume	50–2500 ml
Pressure Support	$0-60 \text{ cmH}_20 \text{ (PEEP+PSV < 100 cmH}_20 \text{)}$
PEEP/CPAP	0-45 cmH ₂ 0
Respiratory Rate	2–120 bpm
Inspiratory Pressure	$5-80 \text{ cmH}_20 \text{ (Peep+Pins} < 100 \text{ cmH}_20 \text{)}$
Inspiratory Pressure Limit	5-80 cmH ₂ 0
Inspiratory Flow	(-) 6–120 L/min (Adult), (-) 6–30 L/min (Paediatric and Neonatal)
Inspiratory Time	0.1-9.9 sec
I:E Ratio	1:9.9–2.5:1
Rise Time (in PCV)	Fast, Medium, Slow (0–2 sec)
Apnea Time	2-60 sec
Apnea Backup Type	PACV, VACV Mode, backup with insp. pressure, Tv
Trend	72 hrs
Enable inhalation trigger (En-sense)	10–80% of Inhaled volume
Exhalation trigger Sensitivity (Ex-trig)	10–80% of peak flow in exhalation
Trigger Sensitivity	Pressure Trigger : OFF, (-)0.5 to (-)20 cmH ₂ 0 Flow Trigger : OFF, 0.5 to 20 lpm
FiO ₂ %	21%-100%
Suitable Body Weight	Volume Mode: 10-150 kg (Paed., Adults), Pressure Mode: 1-150 kg (Neo., Paed., Adults)
Nebulizer (time)	10–180 mins (3 Hours)
Altitude Volume Compensated	(-)100–6560 [m]
Alarms	Airway Pressure / Tidal Volume / Minute Volume / Respiratory Rate / Apnea / O_2 % / Airway Leak / Circuit Open / Battery Indicator / O_2 Supply Pressure Failure / Ventilator In-operative
Power	AC: 100–240V 50–60 Hz, DC : 10.8V, 8.8APower Consumption : 150W
Power Backup	4 Hrs. internal battery backup
International Standard	IEC 60601-2-12:2001 - lung ventilator

Physical Specifications

Dimensions (W x D x H)	250mm x 156mm x 371mm
Weight	7.5 kgs. (main unit with battery)

Environmental Requirements

Storage Temperature	(-) 20°C–50°C
Relative Humidity	10–90%
Operation Temperature	0°C–50°C

Optionals

12.1" TFT Graphical Display Unit	BIPAP Mode
Nebulizer	APRV Mode
Humidifier	Extended Battery (up to 4 hrs)
SpO ₂ Monitoring	Carry Bag
EtCO ₂ Monitoring	Stand



- ,

Meditec International England Limited

8, Pinner View, Harrow, Middx, HA1 4QA, United Kingdom Tel: +44 203 6423808 E-mail : sales@meditecengland.co.uk Website : www.meditecengland.co.uk

#55000315R4 Ver. 19/1 *In view of continued product improvement Meditec International England Limited reserves the right to make changes in the equipment which may affect the information contained in this leaflet without giving any prior notice. Pictures are representative, accessories are optional and may vary at the time of supply.