

ZBS-6E DISINTEGRATION TESTER

The ZBS-6E disintegration tester is specifically designed for disintegration testing of tablets, suppositories, capsules and other solid dosage medications. The main technical indicators complies with USP,BP, JP requirements and specifications about disintegration test.

Features

- Meet Disintegration testing requirements of USP,CHP
- Three baskets may be controlled at the same time, 3 groups test can be carried out simultaneously.
- Standard 2.0mm or 0.425mm diameter stainless steel screen and the necessary block.
- The raising and lowering movement of hanging baskets can start or stop by lifting button
- When test time arrive the set value, the basket will be automatically stopped at the highest location so that basket and beaker are fitted easily.
- Four digital LCD display time
- Set time, Automatic timing , Reset time

Specification

- Timing Range: 0 ~ 24hour
- Predetermined Temperature Range: Room temperature to 45.0°C
Display resolution is 0.1°C
- Accuracy of Temperature Control: $\pm 0.5^{\circ}\text{C}$
- Basket lifting frequency: 30 ~ 32 stroke/min
- Basket Length of Stroke: 55mm \pm 1mm
- Minimal Distance Between Mesh and Beaker Bottom : 25mm \pm 2mm
- Dimension: 420mm x 310mm x 500mm



Tianda Tianfa Technology Co.,Ltd

Address: No.92,Weijin Road,Nankai District,Tianjin,300072,China

Website: <http://www.cntdtf.com> Email: tdtfst@yahoo.cn info@cntdtf.cn

Tel: 86 22-27406167 Fax: 86 22-27401850 mobile: 86 13622069625