



Molecular Sieve Powder 3A

Property	Unit	3A-AG	3A-45	3A-60	3A-D-38	Note
Static H ₂ O Adsorption	%wt	≥25.50	≥25.50	≥23.50	≥25.00	RH75%, 25°C
Exchange Rate	%mol	≥40.00	≥45.00	≥60.00	≥38.00	Flame Photometer
Bulk Density	g/ml	≥0.70	≥0.70	≥0.70	≥0.78	Tapped
Package Moisture	%wt	≤21.00	≤21.00	≤21.00	≤21.00	550°C, 2hr
Screen Residue	%wt	≤0.50	≤0.50	≤0.50	≤0.50	325 Mesh
Cl ⁻	%wt	≤0.10	≤0.10	≤0.10	≤0.10	Colorimetry
PH Value	~	≤11.0	≤11.0	≤11.0	≤11.0	5% solution
D50	um	3.0-5.0	3.0-5.0	3.0-5.0	6.5-9.5	Laser Particle Size Analyzer

Packing:

- 40 kg/bag, and each bag has inner PE bag.
- 500 kg/bag, 600 kg/bag, 1000 kg/bag, and each bag has coating.
- Other packing according to your requirement.