



Molecular Sieve Powder 5A

Property	Unit	5A-75	5A-90	5A-D-83	Note
Static H ₂ O Adsorption	%wt	≥27.00	≥27.00	≥27.00	RH75%, 25°C
Exchange Rate	%mol	≥75.00	≥90.00	≥83.00	Titrimetry
Bulk Density	g/ml	≥0.70	≥0.70	≥0.75	Tapped
Package Moisture	%wt	≤21.00	≤21.00	≤21.00	550°C, 2hr
Screen Residue	%wt	≤0.50	≤0.50	≤0.50	325 Mesh
Cl ⁻	%wt	≤0.10	≤0.10	≤0.10	Colorimetry
PH Value	~	≤11.00	≤11.00	≤11.00	5% solution
D50	um	3.0-5.0	3.0-5.0	6.5-9.5	Laser Particle Size Analyzer

Packing:

- 40 kg/bag, and each bag has inner PE bag.
- 500 kg/bag, 600 kg/bag, 1000 kg/bag, and each bag has coating.
- Other packing according to your requirement.