



Molecular Sieve Powder 13X

Property	Unit	13X	13X-D	Note
Static H ₂ O Adsorption	%wt	≥32.00	≥32.00	RH75%, 25°C
CO ₂ Adsorption	%wt	≥23.00	≥23.00	250mmHg, 25°C
Bulk Density	g/ml	≥0.60	≥0.68	Tapped
Package Moisture	%wt	≤23.00	≤23.00	550°C, 2hr
Screen Residue	%wt	≤0.50	≤0.50	325 Mesh
PH Value	~	≤11.00	≤11.00	5% solution
D50	um	3.0-5.0	5.0-8.0	Laser Particle Size Analyzer

Packing:

- 40 kg/bag, and each bag has inner PE bag.
- 500 kg/bag, 600 kg/bag, 1000 kg/bag, and each bag has coating.
- Other packing according to your requirement.