

FT-2000SE FRIABILITY TESTER

Friability is what reflect an important indicator of the tablet quality. It has a direct impact in the tablet production, packaging and transportation. FT-2000SE friability tester is a specialized equipment which test tablet friability. It is fully functional, easy to operate and its performance are fully in line with the relevant requirements of Chinese Pharmacopoeia and United States Pharmacopoeia.

Features

- According with friability testing provisions of USP,CHP
- Timing and counting function
- Countdown function of timing and counting
- Repeat, Pause and Reset function
- Adjusting speed function
- With wide voltage
- Digital LCD display
- Set time, Automatic timing , Reset time

Specification

- Number of Drum: 2 drums
- Rotational speed Range: 20 ~ 70rpm \pm 1rpm adjustable (Factory preset 25rpm)
- Timing Range: 0 ~ 9hours59minutes59secs alternative countdown
- Counting Range: 0 ~ 9999 times alternative countdown
- Wheel Drum Size: inside diameter: about 286mm, depth: 39mm
- Power Source: (110 ~ 220) \pm 10%V 50Hz ~ 60Hz
- Environment Temperature: 5.0 $^{\circ}$ C ~ 35.0 $^{\circ}$ C
- Relative Humidity: < 80%
- Total Power: 40 W
- Dimension: 420mm x 300mm x 350mm
- Weight: 13kg



Tianda Tianfa Technology Co.,Ltd

Address: No.92,Weijin Road,Nankai District,Tianjin,300072,China

Website: <http://www.cntdtf.com> Email: tdtfs1@yahoo.cn info@cntdtf.cn

Tel: 86 22-27406167 Fax: 86 22-27401850 mobile: 86 13622069626