



Model

Crisp - i

Simple & Captivate





Ceiling Mount

Diagnosing < Surgical


Skill Lab Set up <



Wall Mount



90° straight Binocular Eyepiece Tube





0°-210° Tiltable Binocular Eyepiece Tube



Optics	Full Apochromatic Galilean Type
Eye Piece	10 X Wide Field & High Eye Point (Diopter Adjustment +6 to -6)
Binocular	Straight 90° binocular Optional: 0°210° tilt able binocular tube
Inter Pupillary Adjustamne nt	50mm to 75mm
Magnification	5 Step: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X
Working Distance	f=250mm (Optional 200mm / 300mm / 400mm) Or 200mm – 350mm continuous variable distance manual adjustment
Optical Head Movemnts	0° to 90° axially and 180° side wise
Fine FOcus	Manual Adjustment
Illumination	Homogenous Co-axial shadow free Bright White Spot day light color temperature LED integrated Inbuilt UV filter



Optional: Full HD Integrated Camera

EMC 2014/108/EC EAR Regulation (EU) 2017/745 on Medical Devices Clinical Evaluation report as per EU MDR 2017/745 Technical documentation as per Annex II of EU regulation 2017/745 ISO 13485:2016 Certified facility



Sanma Medineers Vision Private Limited Chennai - 600080, TN, India. sales@sanmas.com www.sanmas.com

Specifications & Design are subject to change without prior notice as a result of further product development ©Samma Medineers Vision Private Limited ALIR ights Reserved