

## Collograft Sterile DRY COLLAGEN SHEET

### ADVANTAGES :

- | Ready to use transparent self Adhesive and Dry membrane that offers complete flexibility .
- | Accelerates Healing rate - minimises dressing.
- | Reduces scarring by depositing oriented and organised fibres
- | Regulates the amount of collagenase expressed by keratinocytes.
- | Reduces wound treatment costs by 40%

Sizes : 5 X 5 cms, 10 X 10 cms, 10 X 25 cms, 15 X 30 cms

### APPLICATIONS :

- | Partial and full thickness wounds
- | Pressure ulcers
- | Venous Ulcers
- | Chronic Vascular ulcer
- | Diabetic ulcers
- | Trauma wounds (abrasions, lacerations, second degree burns, skin tears)
- | Surgical wounds (Donor sites, Grafts, Post moh`s surgery, Post laser surgery, Paediatric and Wound dehiscence)



## Collograft-B

### Dry Collagen Sheet with Non Adherent absorbant Backing

- | Purified TYPE-1 Collagen with Medical Grade PU Foam Backing
- | Flexible and Tacky.
- | Facilitates Healing and Hemostatic
- | PU Foam layer absorbs the exudes and keeps the wound free from contamination.
- | Collagen is absorbed by the tissue in 6 to 10 days



### Applications:

All type of Surface Wounds  
Moderate to Heavily Draining Wounds

Sizes available

5 X 5 cms, 10 X 10 cms, 10 X 20 cms, 15 X 30 cms